



INFORMATIVE NOTE

Name: REFLEXX 36

Date: 01/04/2020 Rev: 017

MANUFACTURER

REFLEXX SpA
Via Passeri, 2 46019 VIADANA (MN)
Tel.: + 39 0375 758891 Fax + 39 0375 464504 info@reflexx.com
www.reflexx.com

Imballo - Pack

100		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
R36 S	8032891631430	S	8032891636435
R36 M	8032891631447	M	8032891636442
R36 L	8032891631454	L	8032891636459
R36 XL	8032891631461	XL	8032891636466

DESCRIPTION: ■ Disposable examination glove not sterile powder free ■ **Material:** Vinyl (PVC) ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** white color ■ **Feature:** smooth surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 1.5 (level G1) for major defects and minor defects

CONFORMITY AND PERFORMANCES

Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)
Reg UE 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 VIRUS (ISO 16604:2004). - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device CLASS I

Reg UE 2017/745 – EN 455 1,2,3 & 4
(CND T01020201 NOT SURGICAL VINYL GLOVES)
- Complies with ASTM-F-1671 (viral penetration test) and ISO 10993-10 regulations (skin irritation and sensitization)



Reg CE 1935/2004 – Reg UE 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length min. (mm)	Width mm (±5)
R36/ S (7)	240	85
R36/ M (8)	240	95
R36/ L (9)	240	106
R36/ XL (10)	240	116


Thickness mm	Palm (+/- 0.02)	0.07
--------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥3.6 N
------------------	------------------	--------

EN ISO 374-1:2016				
Chemical compound :	EN 16523-1:2015 TYPE C			EN 374-4:2013
	COD	Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	2	>30'	52%
Hydrogen peroxide 30%	P	2	>30'	17%
(diethylene glycol mono butyl ether)		2	>30'	79%
Ethylene Glycol Monobutyl Ether(2-butoxyethanol)		2	>30'	73%
Hydrochloric acid 15%		2	>30'	82%
Citric acid 10%		3	>60'	67%
Phosphoric acid 40%		2	>30'	95%
Potassium Hydroxide 15%		3	>60'	41%
Sodium Hydroxide 15%		3	>60'	43%

REG EU 2016/425

EN ISO 374-1:2016/Type C



EN ISO 374-5: 2016



EN 420:2003 + A1:2009

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

Code chemical compound		Protection from
K Sodium hydroxide 40%	CAS:1310-73-2	Protection from use of inorganic bases
P Hydrogen peroxide 30%	CAS: 7722-84-1	Protection from use of peroxides

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product does not contain natural rubber latex ■ The product does not contain nickel ■ In case of transient skin reactions and hypersensitivity reactions, discontinue use and consult a doctor **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 5 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Product Declaration of Conformity on www.reflexx.com