


INFORMATIVE NOTE
Name: REFLEXX 76 soft
Date: 25/10/2022 Rev: 012
MANUFACTURER

 REFLEXX SpA
 Via Passeri, 2 46019 VIADANA (MN)
 Tel.: + 39 0375 758891 Fax + 39 0375 464504 info@reflexx.com
www.reflexx.com

Imballo - Pack

100		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
R76 XS	8032891631874	XS	8032891636879
R76 S	8032891631881	S	8032891636886
R76 M	8032891631898	M	8032891636893
R76 L	8032891631904	L	8032891636909
R76 XL	8032891631911	XL	8032891636916

DESCRIPTION: ■ disposable examination glove powder free, not sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** blue ■ **Feature:** internally chlorinated; textured surface and textured fingertips ■ composition soft and sensitive ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES
Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)

Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg EU 2017/745 – EN 455 1,2,3 & 4


(CND T01020204 NON SURGICAL NITRILE GLOVES)



Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

DIMENSIONS		
Cod./Size	Length (mm)	Width mm (±5)
R76/ XS (6)	230*	75
R76/ S (7)	230*	85
R76/ M (8)	240	95
R76/ L (9)	240	105
R76/ XL (10)	240	115
* Minimum length 240 mm		

Thickness mm	Palm (+/- 0.02)	0.07
EN 455-2:2009/A1	Tensile strength	≥6 N
EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi) pass	Virus protection ISO 16604:2004 pass

REG EU 2016/425

EN 420:2003 + A1:2009

EN ISO 374-1:2016				
Chemical compound:	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	7,1%

RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD	Protection
K Sodium hydroxide 40% CAS:1310-73-2	Protection from use of inorganic bases

RESULTS OF TESTS FOR ADDITIONAL CHEMICAL COMPOUNDS ACCORDING TO EN ISO 374-1: 2016 THE FOLLOWING RESULTS DO NOT GIVE PRESUMPTION OF CONFORMITY:

EN ISO 374-1:2016				
Chemical compound :	COD	EN 16523-1:2015		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Formaldehyde 37%	T	3	>60'	16,02%
Hydrogen peroxide 30%	P	1	>10'	-3,89%

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time

WARNINGS Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In case of transient skin reactions and hypersensitivity reactions, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment

STORAGE AND DISPOSAL Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on www.reflexx.com