



INFORMATIVE NOTE

Name: REFLEXX 77

Date: 17/03/2021 Rev: 07

MANUFACTURER

REFLEXX SpA
Via Passeri, 2 46019 VIADANA (MN)
Tel.: + 39 0375 758891 Fax + 39 0375 464504 info@reflexx.com
www.reflexx.com



Cod./Size	Pack	
	Box	Carton
R77100F/S	8032891633878	8032891638873
R77100F/M	8032891633892	8032891638897
R77100F/L	8032891633922	8032891638927
R77100F/XL	8032891633953	8032891638958
R77100/XS	8032891633465	8032891638460
R77100/S	8032891633472	8032891638477
R77100/M	8032891633489	8032891638484
R77100/L	8032891633496	8032891638491
R77100/XL	8032891633502	8032891638507
R77/XS	8032891633809	8032891638804
R77/S	8032891633823	8032891638828
R77/M	8032891633847	8032891638842
R77/L	8032891633861	8032891638866
R77/XL	8032891633885	8032891638880

DESCRIPTION: ■ disposable glove powder free, not sterile ■ **Material:** nitrile ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** light blue ■ **Feature:** internally chlorinated; textured surface and textured fingertips ■ **Packing:** 100 or 200 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES

Personal Protective Equipment cat III

0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)
Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 VIRUS (ISO 16604:2004). - EN 374-2:2014 - EN 420:2003 + A1 2009

* Medical Device

Reg EU 2017/745
(CND T01020204 NON SURGICAL NITRILE GLOVES)
-* Complies with ASTM-F-1671 (viral penetration test) and ISO 10993-10 regulations (skin irritation and sensitization)



Reg CE 1935/2004 – Reg EU 10/2011 - DM 21/03/73. To such use carefully read the declaration of suitability for food contact (DICA).

* These conformity are not declared on products sold as a food line.

DIMENSIONS				
Cod./Size			Length min(mm)	Width (±5 mm)
/	R77100/ XS (6)	R77/ XS (6)	240	75
R77100F/S (7)	R77100/ S (7)	R77/ S (7)	240	85
R77100F/M (8)	R77100/ M (8)	R77/ M (8)	240	95
R77100F/L (9)	R77100/ L (9)	R77/ L (9)	240	106
R77100F/XL (10)	R77100/ XL (10)	R77/ XL (10)	240	116

Thickness mm	Palm (+/- 0.02)	0.06
---------------------	-----------------	------

EN 455-2:2009/A1	Tensile strength	≥6 N
-------------------------	------------------	------

EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
	pass	pass

EN ISO 374-1:2016				
Chemical compound :	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	7,1%



RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

Chemical compound :		Protection
K Sodium hydroxide 40 %	CAS:1310-73-2	Protection from use of inorganic bases

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ In rare cases where transient hypersensitive skin reactions may occur, discontinue use and consult a doctor ■ The user must inform the manufacturer and the relative authority of any incidents occurring during the use of the protective equipment **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Product Declaration of Conformity on www.reflexx.com