



INFORMATIVE NOTE

Name: REFLEXX 44 black

Date: 18/01/2021 Rev: 07

MANUFACTURER

REFLEXX SpA
Via Passeri, 2 46019 VIADANA (MN)
Tel.: + 39 0375 758891 Fax + 39 0375 464504
info@reflexx.com
www.reflexx.com



Imballo - Pack

100		10	
Taglia Size	Articolo Article	Taglia Size	Cartone Carton
R44 S	80020911732182	S	8772681632127
R44 M	80020911732900	M	8772681632700
R44 L	80020911732915	L	8772681632710
R44 XL	80020911652012	XL	8772681632727

DESCRIPTION: ■ disposable examination glove powder free, not sterile ■ **Material:** natural latex ■ **Shape:** ambidextrous; five fingers, beaded cuff ■ **Colour:** black ■ **Feature:** internally and externally chlorinated; textured surface ■ **Packing:** 100 pieces/box; 10 boxes/carton ■ **ISO 2859:** AQL 1.5 (Level G1) freedom from holes; AQL 2.5 (level G1) major defects; AQL 4 (level G1) minor defects

CONFORMITY AND PERFORMANCES

Personal Protective Equipment cat III0465 ANCI_CIMAC, Via Aguzzafame, 60/b 27029 Vigevano (PV)
Reg EU 2016/425 - EN ISO 374-1:2016 TYPE C - EN ISO 374-5:2016 - EN 374-2:2014 - EN 420:2003 + A1 2009

Medical Device

Reg EU 2017/745 - EN 455 1,2,3 & 4
(CND T010201 NON SURGICAL LATEX GLOVES)

CE
class I
cat III

DIMENSIONS			Thickness mm	Palm (+/- 0.03)	0.12
Cod./Size	Length. min(mm)	Width. (±5 mm)			
R44/ S (7)	240	85			
R44/ M (8)	240	95			
R44/ L (9)	240	106			
R44/ XL (10)	240	116			
			EN 455-2:2009/A1	Tensile strength	≥6 N
			EN ISO 374-5:2016	Protection from microorganisms (bacteria and fungi)	Virus protection ISO 16604:2004
				pass	pass

EN ISO 374-1:2016				
Chemical compound:	COD	EN 16523-1:2015 TYPE C		EN 374-4:2013
		Performance level	Breakthrough time	% degradation
Sodium hydroxide 40%	K	3	>60'	5%



RISKS FROM WHICH THE DPI IS DESIGNED TO PROTECT:

COD	Protection
K Sodium hydroxide 40%	CAS:1310-73-2 Protection from use of inorganic bases

RECOMENDATIONS OF USE Multi-purpose glove to protect the hand and/or the handled product, can also be used for general uses ■ The glove is suitable for contact with pure chemical products for a period of less than the permeation time indicated in the chemical resistance table ■ Results obtained under laboratory conditions may not reflect actual working time **WARNINGS** Before using the product, try it in the real condition of use ■ The glove is a device subject to wear, so frequent replacement is recommended and in all cases where imperfections are observed ■ Wear gloves with dry, clean hands ■ The product contains natural latex which can cause allergic reactions of type I and IV hypersensitivity. In case of reaction, discontinue use and consult a doctor **STORAGE AND DISPOSAL** Store the gloves in their original packaging in a dry place. ■ Do not expose to direct sunlight and heat sources. ■ The expiry date is 3 years from the date of manufacture indicated on the sales package, and is valid for the product properly stored ■ Gloves must be disposed of according to local regulations ■ Declaration of Conformity on www.reflexx.com